

PH (0 - 2.5) BTE 50 BANDELETTES


Référence

0299009

pH 0-2.5 Test Strips

1. Remove one test strip from the bag, holding one end of the strip between your index finger and thumb. Do not touch the test pad with your fingers.
2. Dip the test strip about 1/2" into the solution to be tested for 1-2 seconds.
3. Remove and shake off any excess liquid.
4. Compare immediately to the color chart. Do not allow the test strip to dry before reading against the color chart.

pH range



Disclaimer: Precision Laboratories products are not approved for use as medical devices. Products and product information is for informational purposes only and not intended for use to diagnose, treat, cure or prevent any medical condition. Consult a medical doctor when needed.

Ces spécifications ne peuvent être considérées que comme des informations générales et ne dispenseront pas l'utilisateur de son obligation de s'assurer que le produit est adapté pour l'utilisation prévue. En cas de doute, veuillez nous contacter.

FT 0299009

avril 24